



# BactoSense™ Multi

Automated flow cytometer for easy monitoring of bacteria in multiple water samples



- Autonomous** Processes up to 30 samples automatically, while thermally stabilized
- Safe** Minimized handling of chemicals, all reagents are sealed in a recyclable cartridge
- Accurate** Flow cytometry technology allows precise detection of more than 99.9% of microbial cells
- Fast** Results available within 25 minutes, analyses up to 36 samples in 24h
- Reliable** Self-check routines, factory calibration and low maintenance

- Easy to use** Fully automated sample preparation, batch measurements and cleaning – can be used by anyone
- Cost saving** Reduces the required number of plating tests (HPC) for a minimal total cost of ownership
- Universal** For process monitoring, lab analysis, manual or batch processing – provides TCC, ICC, HNAC/P, and LNAC
- Compact** Built for process and industrial operations
- Integrated** Choice of multiple data interfaces: USB, Ethernet, FTP

## Main applications

Microbiological assessment of your production and processes.

- Allows a complete mapping of your system
- Enables monitoring of multiple sampling points
- Simultaneous disinfection control
- Filtration efficiency
- Synchronized validation of the distribution network

## Industries

- Water treatment & distribution
- Food & beverage
- Laboratories & universities

## Parameters provided

- TCC Total Cell Count
- ICC Intact Cell Count
- HNAC High Nucleic Acid Count
- LNAC Low Nucleic Acid Count
- HNAP High Nucleic Acid Percentage

# Specifications

Measuring principle	Flow cytometry
Light source	Laser diode 488 nm
Optical detection	Fluorescence: 535/43 (FL1), 715 LP (FL2), Side scatter 488/10 (SSC)
Lower size detection limit	0.1 µm
Measuring range	1'000 - 2 Million cells/ml
Detection limit	100 - 5 Million cells/ml
Accuracy	< 5 % relative
Measuring interval	40 minutes
Microbial parameters	TCC/ml, ICC/ml, LNA/ml, HNA/ml, HNAP(%)

Sampling	Autoloading station
Sample volume	1162 µl sampled, 90 µl for analysis
quantity	1 - 30 samples
cooling	min. 4°C, max. room temp. minus 3°C
Requirements	
chlorine concentration	max. 3 mg/l
turbidity	max. 10 FTU
pH range	5 - 12
temperature range	5..40°C
conductivity	max. 100'000 µs/cm at 20°C


Instrument	Factory calibrated
Display	Touchscreen
Data storage	32 GB
Dimensions (WxDxH)	700 × 757 × 373 mm
Weight	35.5 kg
Power supply	100 - 240 VAC, 50/60 Hz, 1.4 A, 2 sockets
Power consumption	20 + 160 W
Ambient temperature	5..30°C
Relative humidity	10 - 90% RH
Cartridge	Hermetically sealed enclosure for reagents, cleaning liquids and waste
Cartridge capacity	Max. 500 measurements, 9 months validity

Interfaces	Digital
Digital interfaces	Sealed USB, Ethernet connection, FTP

Accessories	
TCC Refill	Filling and servicing of cartridge - to measure Total Cell Count of up to 500 samples
ICC Refill	Filling and servicing of cartridge - to measure Intact Cell Count of up to 500 samples
Cleaning kit	Deep cleaning of all internal micro-fluidic components
Validation kit	Easy way to check your instrument after transport or long period out of use
Screw cap vials	10 ml, box of 100 vials with septum
Barcode Reader	Eliminates input errors, reads 1D-, 2D- and QR codes



BactoSense Multi with cooling capabilities



Dotplots showing TCC and ICC



ICC or TCC cartridge



Sample tray containing up to 30 vials



Validation & Cleaning kits



**bNovate Technologies SA**  
 Ch. Dent d'Oche 1A · CH-1024 Ecublens  
 Tel. +41 (0)21 552 14 21  
 info@bnovate.com · www.bnovate.com