BactoSense Quick Guide



This Quick Guide helps you get started with BactoSense, your solution for the automated measurement of the concentration of microorganisms in water.

Before getting started



Follow the user manual instructions

- Read the user manual before installing or operating the instrument.
- Pay particular attention to the safety information in the safety chapter and throughout the manual.
- Make sure you always have access to the user manual.
- Find the latest version of the instrument and accessory user



www.bnovate.com/docs



Verify the operating conditions

- Temperature 5–30 °C (41–86 °F)
- Humidity 10–90 % RH
- Indoor use, max. 2000 m (6600 ft) above sea level



Verify the sample conditions

- In auto mode with the online sampling device at the inlet: flow rate 200– 400 mL/min (0.05–0.10 gal/min), pressure max. 0.5 bar (7.3 psi)
- Temperature 5–40 °C (41–104 °F)
- Turbidity 0–10 FTU, 0–10 FNU
- Conductivity 0–100,000 μS/cm at 20 °C
- Chlorine max. 3 mg/L, iron max.
 200 μg/L, manganese max. 50 μg/L
- pH 5–12

Measuring in auto mode



To measure automatically and in predefined intervals from a connected water source:

- Make sure the correct online sampling device is installed and the sample conditions are met.
- In the *Home menu* of the BactoSense user interface, press *Auto mode* and *Start*.
- On the time series graph view, press *Start*.
- Select the *Online Sampling Device (with Backflush)* protocol and configure the measurement series.
- Follow the instructions on the screen to start measuring.

Measuring in manual mode



- Make sure the manual sampling device is installed correctly.
- In the *Home menu* of the BactoSense user interface, press *Manual mode* and *Start*.
- Select the *Water analysis* protocol and configure the measurement.
- Attach a 5-mL sample tube containing a 2–5 mL water sample to the sampling device.
- Follow the instructions on the screen to start the measurement.

Reviewing results

 In the Home menu of the BactoSense user interface, press Results (Manual) or Results (Auto) and select the desired measurement.

Handling cartridges



The *Cartridge* contains all the reagents needed for measurements and self-cleaning.

• Icons in the top bar of the user interface indicate the filling level of a *Cartridge*.



- Find the expiration date and storage conditions on the cartridge or in the user interface under *System Info*.
- Arrange a *Cartridge* refill once the cartridge level is low or the cartridge is about to expire.
- To install, remove or exchange a Cartridge, press on Maintenance → Cartridge change on the BactoSense user interface and follow the instructions on the screen.

Relocating or storing your BactoSense



Careful decommissioning, transport, storage, and recommissioning is essential. Consult the user manual for this.

Qualification

BactoSense is calibration free.

For operating environments requiring system suitability testing, use the *System Suitability Test Kit* (SST Kit).

For other environments, an instrument qualification using the *Validation kit* is recommended in the following situations:

- After shipping the instrument or other long-distance transport, exposing the instrument to vibration or shock.
- After storing the instrument for more than one month.
- If there are indications that measurement results could be erroneous.

Keeping your BactoSense clean



The Wash Station keeps the manual sampling device clean.

Use the Wash Station:

- Before switching BactoSense off.
- After measuring a highly concentrated sample.
- If the instrument will be idle for more than six hours.
- If the manual sampling device has been newly installed.



The *Cleaning Kit* is used to clear biological and non-biological contaminations.

Use the Cleaning Kit:

- When facing persistent contamination.
- If the needle of the manual sampling device has visible deposits.
- After stopping in an error state with a loaded sample.
- If *High baseline* or *Air bubble detection* warnings frequently appear.

Interpreting user interface icons



The last measurement triggered an alarm.



Maintenance is required. Arrange a maintenance visit.



No cartridge or not initialised correctly. Install the cartridge with on-screen instructions.



Critical error. Resolve the issue and then manually clear the error in the error log.

Contact us

bNovate customer service: <u>support@bnovate.com</u> or your local customer service partner.