

Document number: 40501-01_EN Version: 1 Valid from: BactoSense SW 2.2

Instruction Manual

BactoLink



Connectivity Module

© bNovate Technologies SA, subject to technical changes without notice 10/2023.

This document has been written by bNovate Technologies SA. Copying or modifying the content or giving this document to third parties is permitted only with the expressed consent of bNovate Technologies SA.

bNovate Technologies SA Ch. Dent d'Oche 1A CH-1024 Ecublens Switzerland Tel. +41 21 552 14 21 info@bnovate.com www.bnovate.com

Table of Contents

1 General user information
1.1 Intended Use
1.2 Pictograms
1.3 Overview BactoLink
1.3.1 Prerequisites 2
1.3.2 Power-over-Ethernet (PoE) Injector 3
1.3.3 Antenna 4
2 Installation
2.1 Instrument Connections5
2.2 Connection of Antenna and PoE-Injector
2.3 Mounting BactoLink9
2.4 First installation of BactoLink
2.5 Reconfiguration of BactoHub14
3 Remote Connections15
3.1 Web User Interface (WebUI)
3.2 Graphical User Interface Relay (GUI relay)
4 Troubleshooting

General user information

1.1 Intended Use

BactoLink is intended to function alongside BactoSense and BactoSense Multi by establishing a remote connection between the instrument and your PC. This connection is initiated through the BactoHub portal using HTTPS and your web browser. The BactoHub portal serves as the access point for BactoLink and the associated BactoSense devices. It's important to note that the connection to your instrument remains active for 30 minutes and must be re-established after this period.

1.2 Pictograms

All pictograms used in this document are explained below:



Additional information about the current topic.



Practical procedures when working with the BactoSense.



The screenshot is an example and may differ from the current device.

1.3 Overview BactoLink

BactoLink consists of an antenna, an ethernet cable and a Power-over-Ethernet (PoE) injector. The PoE injector is connected to the Power Supply of the BactoSense (Multi).

1.3.1 Prerequisites

Before using BactoLink, ensure the following prerequisites are met: The instrument has at least Software Version 2.2 or above. DHCP and HTTPS must be activated within the Instrument settings. Optionally, if GUI relay is required activate VNC in the instrument settings. Ensure that you have access to an account at hub.bnovate.com. If a new account is required you can request one at support@bnovate.com.

1.3.2 Power-over-Ethernet (PoE) Injector

The PoE-Injector delivers power to the antenna through an Ethernet connection. Both the Ethernet cable connection and the power connection to the BactoSense (Multi) are permanently integrated into this device. It includes a port for the power cable provided with the BactoSense (Multi) Power Supply unit. The final connection involves linking the Ethernet to the antenna. For comprehensive installation guidelines, please consult the "Installation" chapter.



Figure 1-1Dimensions of the PoE injector. PoE injector height is 65.5 mm from feet to the top.

1.3.2.1 Mechanical Specifications

Table 1-1 Mechanical Specifications of the PoE injector

Materials	Polycarbonate (PC)	
Connector	RJ45, Power Connector	
Protection	IP67	
Dimensions	165 x 165 x 54 mm (LxHxW)	
Operating Temperature	5°C to 30°C	

1.3.3 Antenna

The antenna is connected to the PoE Injector via an Ethernet cable. The antenna houses the SIM card for mobile connections (4G/LTE, 3G, 2G) to the BactoSense (Multi). Information on the mountain bracket can be found in chapter "Mounting BactoLink".



Figure 1-2Perspectives and dimensions of the antenna

1.3.3.1 Mechanical Specifications

Table 1-2 Mechanical Specifications of the antenna

Materials	ABS, PC, FR4
Connector	RJ45
Protection	IP67
Dimensions	165 x 165 x 54 mm (LxHxW)
Weight	0.8 kg
Operating Temperature	-40°C to 80°C

2 Installation

2.1 Instrument Connections



1	Ethernet connection between PoE-Injector and RJ45 port of BactoSense (Multi)
2	Power connection between PoE-Injector and power port of BactoSense (Multi)
3 Ethernet connection between PoE-Injector and antenna	
4	Power supply to PoE-Injector. Note: The power supply is delivered with BactoSense, not with BactoLink.

2.2 Connection of Antenna and PoE-Injector

	WORK STEP	ADDITIONAL INFO / IMAGE
1.	Open the bag that contains the ethernet cable	
2.	Open the bag that holds the PoE injector. Inside, there should be a second bag containing a connector. Proceed to open this second bag. It should contain the following parts: • Cap • Sealing ring with opening • Spacer • Connector to PoE Injector Line them up on the ethernet cable like shown to the right.	Cap Sealing Ring with Opening Depring Ethernet Cable
3.	 Guide the RJ45 connector into the connector of the PoE injector. Both connectors allow for only one orientation. Push the RJ45 connector as far to the front as possible. Push the spacer and sealing ring into the connector to PoE injector. Screw the Cap on the Connector, gently pushing all the parts on the ethernet cable together. Screw until hand-tight. The resulting fitting is a bayonet fitting to the PoE injector 	
4.	Take the other side of the ethernet cable and align it in front of you.	
5.	In the package that contains the antenna there should be the following parts: • Antenna Connector • O-ring	

	WORK STEP	ADDITIONAL INFO / IMAGE
6.	Unscrew the connector and push the sealing ring out. The curved rod from chapter "Mounting BactoLink" may be used to do this, if this proves hard to do.	
7.	Align the parts of the antenna connector on the ethernet cable like shown on the picture on the right. Note: The critical sequence involves connecting the RJ45 connector to the antenna first, securing the connector to the antenna, and finally, mounting the cap to prevent the cable from being locked beneath the seal.	Cap Antenna Connector to Antenna Connector Con
8.	Push the RJ45 connector lightly into the connector to the antenna. See that the little stub on the RJ45 is within the antenna.	
9.	Push the sealing ring into the connector to the antenna.	

	WORK STEP	ADDITIONAL INFO / IMAGE
10.	Put the O-ring on the connector to the antenna and connect the RJ45 connector.	
11.	Screw the connector to the antenna in the antenna.	
12.	Screw the cap on the connector to the antenna. The final assembly should look like the picture to the right.	

2.3 Mounting BactoLink

BactoLink includes a mounting set for connecting the antenna to cylinder or square-shaped objects, pipes or rods.

\sim		~
(4)
L		1
L		1
U	0	J
\sim		~

	WORK STEP	ADDITIONAL INFO / IMAGE
1.	Open the package containing the antenna for the BactoLink. Inside, you will find a mounting bracket and the necessary components for assembly.	
2.	 Take out the mounting bracket. It consists of three essential parts: A. Mounting Plate: This is a flat piece of metal designed to be attached to the BactoLink. B. Curved Rod: A bent rod holds the BactoLink in place. C. Spacer: A small component maintaining the distance between the curved rod and the mounting plate. 	
3.	Locate the two screws with washers provided in the antenna package. These will be used to secure the mounting plate to the BactoLink.	000
4.	Position the mounting plate on your BactoLink. Insert the screws through the holes in the mounting plate and then into the corresponding holes on the BactoLink. Ensure that the washers are placed between the screw heads and the mounting plate. Using an Allen key (5 mm), tighten the screws until the mounting plate is securely attached to the BactoLink.	Threads for mounting bracket

	WORK STEP	ADDITIONAL INFO / IMAGE
5.	Locate the two flanged nuts provided in the antenna package. These will be used to secure the spacer and curved rod to the mounting plate.	
6.	Now, the curved rod and spacer will be attached to the mounting plate. Position the curved rod and spacer onto the mounting plate in the desired orientation.	Mounting Plate Screws Mounting Plate Screws Flanged Spacer Nuts Spacer Curved Rod A - From Right B - From Left
7.	Insert two flanged nuts onto the threaded end of the curved rod. Hand tighten or use an Allen wrench (13 mm) until no movement of the flanged nuts is possible, and the BactoLink is securely attached.	

2.4 First installation of BactoLink

_	_
(
U	J

	WORK STEP	ADDITIONAL INFO / IMAGE
1.	Check whether the instrument is running on Software Version 2.2 or higher and upgrade if necessary. If you need a software upgrade contact support@bnovate.com.	
2.	To connect to the BactoHub, scan the 2D barcode on the backside of the antenna (picture to the right is exemplary) and enter your BactoHub credentials, composed of the username and password you chose earlier for BactoHub. If you do not have a BactoHub account, contact support@bnovate.com.	Threads for mounting bracket
3.	Power off the instrument and connect the devices like described in Chapters 2.1and 2.2	<section-header><section-header><section-header><complex-block><image/><image/><image/><image/><image/></complex-block></section-header></section-header></section-header>
4.	Power on the device.	
5.	Ensure DHCP is activated. If not activate it. (Home Menu > System Settings > Network) Do not power cycle the instrument immediately.	

	WORK STEP	ADDITIONAL INFO / IMAGE
6.	Ensure HTTPS for access to WebUI is activated and assigned to port 443. If not, activate and reboot BactoSense. (Home Menu > System Settings > System Services) MTTPS is necessary for BactoLink to function.	
	Power cycle the instrument or perform optionally step 7.	
7.	If GUI relay (remote desktop) is needed, ensure VNC is activated and assigned to port 5900.	Connect BactoLink
	(Home Menu > System Settings > System Services)	Manage system service. Detailstand?************************************
	Be aware that with GUI relay your actions on the screen are visible to others, either standing in front of the instrument or remotely connected. Confidential information might be shared.	
	Power cycle the instrument.	

	WORK STEP	ADDITIONAL INFO / IMAGE
8.	Configure BactoLink via the BactoLink serial number and a user-friendly name. Click on "Connect".	
9.	Wait until the instrument is connected, then click on "Finish".	<page-header> but the set of the set</page-header>
10.	The instrument with BactoLink will appear in BactoHub in your instrument list as shown on the right. Two connection possibilities are given: WebUI and GUI relay. Once the connection is established, the window will open in a separate tab. The connection will be valid for 30 min. If this time frame is exceeded, a new connection will need to be made.	<page-header><text><text><section-header><section-header><section-header><section-header><section-header><section-header><section-header><section-header><image/><section-header></section-header></section-header></section-header></section-header></section-header></section-header></section-header></section-header></section-header></text></text></page-header>

2.5 Reconfiguration of BactoHub

	WORK STEP	ADDIT	IONAL INFO / IMAGE	
1.	Reconfigure BactoLink by clicking on the button "Initialize BactoLink" below the "My Instruments" header. Enter the necessary data: Instrument s/n and user friendly name. Press continue.	€ Microsoft Ourstantion My Bochannets Wy Certrisign Relif Research Wy Asthematic	With the state st	Sara Faceli
2.	Finish the steps 2-9 from chapter 2.4. You can now access your devices through WebUI and GUI relay.			

3 Remote Connections

There are two connection possibilities: WebUI and the remote desktop of the GUI relay, which are described in the following chapters.

3.1 Web User Interface (WebUI)



The Web User Interface allows access to data on BactoSense (Multi) via the RJ45 port. Clicking on the WebUI button in the BactoHub establishes the connection via BactoLink. If you reboot BactoSense through the Web User Interface, you need to reload your web page as well.



Password	Enter the credentials when prompted (Username - Basic, Advanced,
	Admin, or service - and the respective chosen Bactosense device
	password).



Dashboard	Provides an overview of the latest measurement results, instrument
	status and Warnings if a Cartridge is expired.



Dotplots	Showcasing a video using recent measurements. It's important to note that loading may require some time, depending on the connection. Therefore, it's advisable to initially set the animation speed to 1. Once the images are fully loaded, higher animation speeds can be used. Additionally, this setup should offer download options in both GIF and
	FCS formats.

PLOT															
															W & 4
					In	strument nam	e:							10 M T 11 II	20 W
14 1	ter ter for ty all												- 100	- HNAP 200	
	(a) (10) (00) (11) (40	243 011	×											had but	
1004															609
3															501
20		35 456321	HNAP							_					40
8 10k															30
0,															20
1000															10
	Oct 2	042 2021 21:19	Oct 4	Oct 6	Oct 8		oct 10	00	tt 12		Oct 14	Oct 14	E Contraction of the second se	Oct 18	
DOWNLOAD RES	SULTS			_	_						-				
DOWNLOAD RES Download all data in Cl • Auto moder dat • Manual moder dat • Instrument valid	SULTS CVr or Encel format: in cvr pdf siter cvr pdf darlee: xiter cvr pdf														
DOWNLOAD RES Deveload al data in C A date mode the Manual mode it Instrument uald	SULTS CSV or Excel format: a fox pdf altar cov pdf darlien: alsa cov pdf SULTS			_	_										
DOWNLOAD RES Develoed al data in Ci A data incide da Manual mode da Instrument uaid	SULTS SV or Excel formal: is (cv) pdf stat (cv) pdf dation: slot (cv) pdf SULTS 21 - 20 October 2023 ~	Download 200													
DOWNLOAD RES Download al data in Co Adda mader data Manual moder da Instaument adda AUTO MODE RES Conter Conter Date	SULTS CSV or Eucel format: (a or) pdf stat (av) pdf danion: ska; (av) pdf SULTS 21 - 20 October 2023 ~ Sampling date ~	Download 208			~ LNAC	< HNAC	 нкар 	× TCC ×	KC - DCC	v	v				
DOWNLOAD RES Devinited all data in C 4 Ada node in 9 Instantineder 9 Instantineder 14 Instantinet auf 14 AUTO MODE RES 14 Zehnung 2022 Date v	SULTS EXV or Bacel Kornut: In care jed skel (care jed) subset (care jed) SULTS 21 - 30 Cataber 2023 ~ Sempling date ~	Downiced 20#			- UNAC	~ HRAC	 HNAP 	× TCC ×	NCC ~ DCC	*	v				
DOWNLOAD RES Develops 34 drist in C 4 Adds mode find 9 Ministration 1 Adds mode 9 Instrument Laid AUTO MODE RES 2012-19-20 3419(10	SULTS EXY or Back Yomat: In can jed Hall (ox) jed danien: Xox (cx) pot SULTS 21 - 30 October 2023 ~ 21 - 30 October 2023 ~ 21 - 30 October 2023 ~ 21 - 20 October 2023 ~ 21 - 20 October 2023 ~ 2023-13-20 Yaža Ma	Download 20#	_	_	~ LNAC 6130	~ HNAC 99440	· HNAP 62	× TCC ×	NCC ~ BCC 161280	•	v ing (cs.				
DOWNLOAD RES Deveload all data in C Adda make india dia data Manada make it Manada make it Matanama data Auto MODE RES 2012 Pelosary 2021 Date 2023-19-20 14 199 10 2023-19-20 11 19 10	SULTS CSV or Suei format: a [cw] pdf skal [cw] pdf skal [cw] pdf SULTS SUL	Downiced 7/# Name			< ENAC 6130 71740	 HNAC 99440 116400 	* HRAP 42 60	× TCC ×	NCC ~ DCC 161280 191140	*	× ing (rs. ing) (s.				
DOWNLCAD 855 Download 34 data av 4 data av	SULTS CX or Seel Konst: a (cs) [cs] stel (cs) (cs) stel (cs) (cs) SULTS 21 - 30 October 2023 ~ Sampling date ~ 2033-19-20 152.54.54 2023-19-20 052.64.54	Download 7/8 Name			 ENAC E110 E110 T783 H110 	 HNAC 99400 19600 19700 	* HHAP 62 60 27	× TCC ×	NCC DCC 191200 191400 223200	•	v ing (to: ing) fo: ing (to:				
DOWNLOAD RES Develop all data olic di Auto and colic di Hintare data olic di Hintare data olic di Hintare data olic di Elebargo di Elebargo data olicitato Elebargo di Elebargo	SULTS Convertisent format: a loss just what i toxi just distributives what i cov just status Sutternative status SULTS Sampling date V 2023-19-00 Head At 4 2023-19-20 Head At 4 2023-19-20 Head At 4 2023-19-20 Head At 4 2023-19-20 Head At 4 Convertisent At 4 2023-19-20 Head At 4 2023-19-20 Head At 4 2023-19-20 Head At 4	Deveload 20 ⁹ Name			 IMAC 4430 7743 4433 1739 	 HNAC 99400 115400 80100 80100 	 HNAP 42 42 42 43 43 44 45 46 46 47 47 48 49 49 49 49 49 49 49 40 <li< td=""><td>× TCC ×</td><td>NCC DCC 191200 191400 223280 2337410</td><td>*</td><td>v Ing (s. ing) (s. ing (s. ing) (s. ing) (s.</td><td></td><td></td><td></td><td></td></li<>	× TCC ×	NCC DCC 191200 191400 223280 2337410	*	v Ing (s. ing) (s. ing (s. ing) (s. ing) (s.				
DOWNLOAD RES Download al data of Anal mode al A	SULTS SUP of Deal format: a (ov) (pdf datalen r/kr (1 ov) (pdf data v) SUP of Deal r/kr (1 ov)	Download 20F			 UNAC 4400 7743 44135 11789 11889 	 NNAC 95400 115400 80150 80150 80170 	 HNAP 482 462 462 37 37 	× TCC ×	NC - DCC 161280 223280 237410 215560	×	v ng ta ng ta ng ta ng ta				
DOWNLOAD RES Deveload all data C Ada made via Monad mode a M	SULTS Chr of boot format: a (ce) [aff shat (ce) [rdf SULTS Sumpling data 2023-19-20 Classes V 2023-19-20 Table 14 2023-19-20 Table 14 2023-19-20 Data 14	Download 20			* UNAC 6030 12784 14130 13890 13880 13880	 HRAC 99400 115400 80300 80300 80310 	• HBA5 42 40 37 37 37 37	× TCC ×	NCC - PCC 141280 194140 223380 223380 223380 22340	*	V Reg (to: reg (to: reg (to: reg (to: reg (to: reg (to:				
DOWNLOAD RES Dreveload at data in C Adda mathin water and a construction of the data in C Adda mathine water andata in C Ad	SULTS Size to boll format: and pipel and pipel and pipel subset. All configures SULTS Sults	Durefield 29			 UMAC 4100 4103 4103 4103 1028 10282 	 HNAC HNAC 9440 116400 8070 8020 80270 8020 8020 	 HBAP HBAP 42 60 37 37 37 	× TCC ×	NCC ~ DCC 161280 151280 2217410 2217410 2215840 2215440 2215410		* * * * * * * * * *				

Results	Shows a plot of measurements over time and provides download options. Below the plot download areas are located for auto mode, manual mode and instrument validation. Additionally, the results of auto mode & manual mode are listed at the bottom of the page. In this overview you can see the results in table format, download cumulative ZIP files and the IMG and FCS files from each measurement. For BactoSense Multi, only the Manual Mode is displayed and no plot of
	BactoSense Multi, only the Manual Mode is displayed and no plot of the measurements over time will be shown.

😑 Statuc Idle.						Ingish Admin
Darbourt	PROTOCOL LOG					
III Results	Date	~ Type	Protocol	< Name	- Success	
III Protocol log	2023-10-20 14:10:10	Auto	analysis_water_		~	1
A tracks	2022-10-20 11:19:10	Acto	analysis, seater,		~	
	2023-10-20 08.19.11	Asto .	analysis_water		~	
Ə system info	2023-10-20 05:19:11	A480	analysis, water,		~	
C. Settings	2023-10-20 02:19:11	Auto	analysis_water		v	
	2023-10-19.22.19:10	Auto	analysis_water_		~	
	2022-10-19.20:19:11	ALED	analysis_water_		*	
	2023-10-19 17:19.10	Auto	analysis_water_		~	
	2023-10-19 14:19:10	Auto	analysis_water_		~	
	2022-10-19 11:19:10	Auto	aralysis_state_		~	
	3023-10-19.08.99.11	Zi den	analosis water		J	•

Protocol Log	Offers an overview of run protocols and error messages. A CSV file
	containing the respective data can be downloaded.

 nois off-strat Cov Line. 		
CONCE LON		
Dute	 Severity < Name 	~ R
1		
2023-09-29 30:29:22	CO1 - Errors cleaned automatically	
2023-09-15 22:02:0	0	
2023-07-25 54:32:04	COI - Errors cleared automatically	
2023-07-02 23:23:19	0	
2022-08-24 13:53:48	C01 - Errors deared automatically	
2022-00-13 14 27:10	0	
2021-10-29 13:58:16	CO1 / Errors cleared automatically	
2021-10-22.00:36:38	o	
2021-03-24 14:27:38	C02 - Errors cleared manually	
2021-03-24 34.8223	0	
WARNINGS		
WARNINGS Dute	× Name	
WARNINGS Date	- Nor	
WARNINGS Date 2003-00-13 201 01:00 2003-00-13 201 01:00	- Nana	
WARNINGS Data Data Data Data Data Data Data Dat	- None	
944ANING5 Den 2003-de 17 2016-0 2003-de 17 2016-0 2013-de 17 2016-0	* Nor	
84ANN055 Per 2009-0412000 2009-0412000 2009-0412000 2009-0412000 2009-0412000	- Nor	
9424091655 Dea 24209-074030 24209-074030 24209-074030 24209-074030 24209-074030 24209-074030 24209-074030 24209-074030 24209-07403 24209-0740 242000 24200 242000 24200 24200 2420	• Nor 11 27 54 20 20	
UAANNOS De 2024 04 13 7891 2024 04 13 7891 2024 04 13 7891 2024 04 13 7891 2024 04 13 7891 2024 04 13 7891 2024 04 13 7891 2024 04 13 7893	- Nori	
94.509/955 Des 2409 49 (1999) 2409 49 (1999	• Nor IV 7 4.4 4.5 4.5 4.5 4.5 4.5 4.5 4.5 4.5 4.5	
BANKHOS Des 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008) 2006-041 (2008)	- Nex	
94.000.965 De 2020-961.958.90 2020-97.958.9	• Non 	
UANNES De 20년 4년 120년 18 20년 4년 1301 18 20년 4년 1301 18	- Net H	

Error Log	Lists all errors and warnings, which can be downloaded as XLSX, CSV
	or PDF.

CARTRIDGE		NETWORK		SYSTEM INFO	
		and a second second			
Cartridge type	KC-A	MAC	00:13:05:3cfb:f7	Device name:	
Cartridge level	42.%	. m	192.168.2.116	Instrument serial number:	
Cartridge serial number		websecket	Connected	Software version	
Cartridge fill number.	BND-ICC			Remaining disk space:	2529.08
Expiry date:	Jun 21, 2024			Erough spice for approx.	0,595 measurements
ourours		TRANSCOMMENT		SERVICE	
and the second se					
Output	Level	Enclosure temperature	21.7 °C	Nest service due:	Nov 23, 2023
Analog 1	2 mA	Enclosure humidity:	24.7 %8H	Value lifetime remaining:	80.16
Analog 1	13.88 mA	Rapid heating on:	0	Plunger lifetime remaining:	45 %
Digital 1	15	Standby heating ox.	1.		
Digital 2	10				
Digital 2	1				
Digital 4	1				
SCREENSHOT					
and the second se					
Get screenshot					

System Info Provides general system information and allows taking screenshots.

🚍 Status: Idle.			English	Admin
🏟 Dashboard				
O Dotplots	NOTOMATIC MODE			
I Results	Sample name	Manual Sampling Device		
🗮 Protocol log	Measurement interval [minutes]	120		
A Error log				
पुर System info		Save		
📽 Settings 🔹				
	INSTRUMENT CONFIGURATION			
	Instrument name	BactoSwin		
		Save		

Settings	Allows changing sample names, measurement intervals, and the		
	instrument name. Not available for BactoSense Multi.		

3.2 Graphical User Interface Relay (GUI relay)

The remote desktop of the GUI relay allows you to access the device as if you were physically present.



Figure 3-1 Exemplary picture of the remote desktop of the GUI relay



Be aware that your actions on the screen are visible to others, either standing in front of the instrument or remotely connected. Confidential information might be shared.

4 Troubleshooting

If BactoSense (Multi) was turned off during installation, BactoLink should work immediately. If issues persist:

- 1. BactoLink must be fully powered on before BactoSense is powered on to establish the connection from BactoSense to BactoLink. Allow at least 3 minutes for the instrument to establish a network connection.
- 2. Check mobile network connectivity at your location using your cell phone.
- **3.** If no connection can be made, disconnect and reconnect the antenna from the PoE injector to power cycle the antenna.
- 4. Inspect all connections on the PoE injector and connected instruments.
- 5. If issues persist, power cycle BactoSense.

bNovate Technologies SA Ch. Dent d'Oche 1A CH-1024 Ecublens Switzerland

Tel. +41 (0)21 552 14 21 info@bnovate.com www.bnovate.com